

Application Brief

Medical Wristbands



Industry
Healthcare

Products
INband Wristband
PF8

Region
North America

Typical Applications
Patient identification

Benefits

- Positive patient identification that helps reduce medical errors
- Ease of application
- High image durability - designed to survive typical exposures in healthcare environment
- Lightweight, flexible construction that ensures patient comfort

Industry Need

Effective patient care is predicated on accuracy, but human error is exceptionally difficult to eliminate in the inherently high-touch medical environment. To address this, barcodes and automatic ID technologies are being used increasingly across a wide range of medical practices. Patient wristbands are one common application area. Intermec patient wristbands ensure positive identification for integration with verification systems for medication, samples, and surgical procedures.

The five rights of medication safety (right patient, right drug, right time, right dose, right route) provide an important framework for verification of pharmaceutical administration, but commonly rely on manual human checks. Integration of positive patient ID with barcode wristbands provides a highly reliable adjunct to eliminate human error. With increasing implementation of patient electronic medical record (EMR), scanning the wristband can trigger caregiver notification of allergies or other special requirements.

The Joint Commission on Accreditation of Healthcare Organizations (JCAHO) has made improving the accuracy of patient identification a high priority, and it remains the number one goal in the pre-release version of the 2010 National Patient Safety Goals (NPSGs) for hospitals. Recognizing the accuracy enhancements of auto ID technologies, bar code verification is specifically cited as an acceptable alternative to resource-intensive two-person verification required for administration of blood products.

Reference Links

JCAHO 2010 National Patient Safety Goals
<http://www.jointcommission.org/PatientSafety/NationalPatientSafetyGoals/>

Intermec Solution

Intermec INband wristbands provide the positive identification hospitals need for maximum patient safety. Excellent print quality ensures both standard and 2D bar codes will be scanned accurately. Durable topcoating provides high resistance to soaps, hand sanitizers, moisture, and other conditions commonly encountered in the hospital environment. Flexible and lightweight, the Intermec band is easy to apply, minimizing the curl that may be experienced with other products. High tensile strength film and a high-shear adhesive closure securely affix the band in place for the duration of the patient stay. A robust tamper-evident closure feature protects the integrity of the patient identification system.

Direct thermal technology delivers mission-critical printing system reliability while providing straightforward compliance with HIPAA regulations for security of patient information. Simple, intuitive printer loading and one-at-a-time printing allows bands to be generated with minimal waste.

Products

INband Wristband

Printer Recommendation

PF8 Desktop Printer [Profile](#)

- *Easy to Deploy* – Features internal Ethernet option and robust configuration tool
- *Easy to Install* – Uses common media and printer languages
- *Easy to Use* – Drop-in media loading

The Intermec Advantage

With safety and quality of care depending on accurate patient identification, wristband system components must work together flawlessly. Intermec printers and media are optimized to deliver superior performance when used together. Our rigorous testing and co-engineering ensures consistently high print quality, proven performance in demanding real-world environments, and maximum print head service life for reduced downtime.

Intermec printers and media are complemented by our complete line of integrated solutions and services, including rugged mobile computers, RFID systems, and advanced scanning technology. Intermec products and services are used by customers worldwide to improve the productivity, quality and responsiveness of their business operations.

For more information

Contact Intermec for more information on positive patient identification with INband wristbands.



North America

Corporate Headquarters

6001 36th Avenue West
Everett, Washington 98203
Phone: (425) 348-2600
Fax: (425) 355-9551

North Latin America

Headquarters Office

Mexico
Phone: +52 55 52-41-48-00
Fax: +52 55 52-11-81-21

South Latin America

Headquarters Office

Brazil
Phone: +55 11 5502.6770
Fax: +55 11 5502.6780

Europe/Middle East & Africa

Headquarters Office

Reading, United Kingdom
Phone: +44 118 923 0800
Fax: +44 118 923 0801

Asia Pacific

Headquarters Office

Singapore
Phone: +65 6303 2100
Fax: +65 6303 2199

Internet

www.intermec.com
Worldwide Locations:
www.intermec.com/locations

Sales

Toll Free NA: (800) 934-3163
Toll in NA: (425) 348-2726
Freephone ROW:
00 800 4488 8844
Toll ROW: +44 134 435 0296

OEM Sales

Phone: (425) 348-2762

Media Sales

Phone: (513) 874-5882

Customer Service and Support

Toll Free NA: (800) 755-5505
Toll in NA: (425) 356-1799



Copyright © 2010 Intermec Technologies Corporation. All rights reserved. Intermec is a registered trademark of Intermec Technologies Corporation. All other trademarks are the property of their respective owners. Printed in the U.S.A. 612047-01A 05/10

In a continuing effort to improve our products, Intermec Technologies Corporation reserves the right to change specifications and features without prior notice.